

MENGHINI PUNCTURE NEEDLE

Ø 1.2 x 70mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-03

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Biopsy Needles |
| Aesculap | SF908 |
| KLS Martin | 18-700-20 |
| Medicon | 09.75.20 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Biopsy Needles |
| Aesculap | SF908 |
| KLS Martin | 18-700-20 |
| Medicon | 09.75.20 |