

MENGHINI PUNCTURE NEEDLE

Ø 1.2 x 100mm; Luer-Lock connection

Minimum order quantity: 10 pieces

SKU: 086-115-04

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF909
KLS Martin	18-700-25
Medicon	09.75.25

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Biopsy Needles
Aesculap	SF909
KLS Martin	18-700-25
Medicon	09.75.25