

## MENGHINI PUNCTURE NEEDLE

Ø 1.4 x 100mm; Luer-Lock connection

Minimum order quantity: 10 pieces

**SKU:** 086-115-06

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Biopsy Needles</a>
<b>Aesculap</b>	SF913
<b>KLS Martin</b>	18-700-35
<b>Medicon</b>	09.75.35

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Biopsy Needles</a>
<b>Aesculap</b>	SF913
<b>KLS Martin</b>	18-700-35
<b>Medicon</b>	09.75.35