

## METAL JARS

Ø 100x200mm



**SKU:** 250-142-09

**Categories:** [Open Surgery](#), [Various Items](#)

### ZUSÄTZLICHE INFORMATIONEN

|                   |                       |
|-------------------|-----------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Bowls</a> |
| <b>KLS Martin</b> | 56-360-99             |
| <b>Medicon</b>    | 47.40.12              |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                       |
|-------------------|-----------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Bowls</a> |
| <b>KLS Martin</b> | 56-360-99             |
| <b>Medicon</b>    | 47.40.12              |