

MICRO-ADSON DRESSING FORCEPS

straight; serrated

15.0cm (6"); width of jaw 1.2 +/-0.1


SKU: 020-076-15

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD223
KLS Martin	12-168-15
AUTHORS	ADSON
Medicon	06.22.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD223
KLS Martin	12-168-15

AUTHORS

[ADSON](#)

Medicon

06.22.15

SPECIFICATIONS

[serrated](#), [straight](#)