

## MICRO-ADSON DRESSING FORCEPS

straight; serrated

15.0cm (6"); width of jaw 1.2 +/-0.1


**SKU:** 020-076-15

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>  |
| <b>TYPE</b>           | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>Aesculap</b>       | BD223   |
| <b>KLS Martin</b>     | 12-168-15   |
| <b>AUTHORS</b>        | <a href="#">ADSON</a>   |
| <b>Medicon</b>        | 06.22.15  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a>               |



### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>Aesculap</b>   | BD223   |
| <b>KLS Martin</b> | 12-168-15   |

## AUTHORS

[ADSON](#)

## Medicon

06.22.15

## SPECIFICATIONS

[serrated](#), [straight](#)