

## MICRO-ADSON DRESSING FORCEPS

straight; serrated

12.0cm (4.8"); width of jaw 0.9 +/-0.1


**SKU:** 020-076-12

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD220
<b>KLS Martin</b>	12-404-12
<b>AUTHORS</b>	<a href="#">ADSON</a> , <a href="#">MICRO</a>
<b>Medicon</b>	06.22.12
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>Aesculap</b>	BD220
<b>KLS Martin</b>	12-404-12

## AUTHORS

[ADSON](#), [MICRO](#)

## Medicon

06.22.12

## SPECIFICATIONS

[serrated](#), [straight](#)