

MICRO DRESSING FORCEPS

straight; serrated; tips 0.8mm

15.5cm (6.2"); width of jaw 0.8mm

SKU: 020-096-08

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
KLS Martin	12-184-08
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ear and Nose Dressing Forceps
KLS Martin	12-184-08
SPECIFICATIONS	serrated , straight