

MICRO DRESSING FORCEPS

straight; serrated; tips 0.8mm

15.5cm (6.2"); width of jaw 0.8mm

SKU: 020-096-08

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Ear and Nose Dressing Forceps |
| KLS Martin | 12-184-08 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Ear and Nose Dressing Forceps |
| KLS Martin | 12-184-08 |
| SPECIFICATIONS | serrated , straight |