

## MICRO DRESSING FORCEPS

straight; serrated; tips 0.6mm

15.5cm (6.2"); width of jaw 0.6mm

**SKU:** 020-096-06

**Categories:** [Forceps](#), [Open Surgery](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear and Nose Dressing Forceps</a>
<b>KLS Martin</b>	12-184-06
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ear and Nose Dressing Forceps</a>
<b>KLS Martin</b>	12-184-06
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>