

## MICRO EAR CUP FORCEPS

straight; 1.8 x 0.9mm

shaft 8.0cm (3.2")



**SKU:** 111-344-08

**Categories:** [Open Surgery](#), [Otology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro Ear Forceps</a>
<b>Aesculap</b>	OG345
<b>KLS Martin</b>	36-003-01
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro Ear Forceps</a>
<b>Aesculap</b>	OG345
<b>KLS Martin</b>	36-003-01
<b>SPECIFICATIONS</b>	<a href="#">straight</a>