

MICRO LARYNGEAL SCISSORS

straight; Ø 2.0mm tube shaft

24.0cm (9.5")



SKU: 131-603-00

Categories: [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Endolaryngeal Micro-Surgery |
| Aesculap | OP310 |
| Medicon | 73.30.40 |
| SPECIFICATIONS | straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Endolaryngeal Micro-Surgery |
| Aesculap | OP310 |
| Medicon | 73.30.40 |
| SPECIFICATIONS | straight |