

## MICRO NEEDLEHOLDER; ROUND HANDLE

straight; 0.3mm; with catch

15.0cm (6"); diamond coated

**SKU:** 096-002-15

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>KLS Martin</b>	20-000-01
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">with catch</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>KLS Martin</b>	20-000-01
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">with catch</a>