

## MICRO NEEDLEHOLDER; ROUND HANDLE

straight; 0.4mm; w/o catch

18.0cm (7.2"); diamond coated

**SKU:** 096-100-18

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>KLS Martin</b>	20-004-18
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">without catch</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Micro-Needle Holders</a>
<b>KLS Martin</b>	20-004-18
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">without catch</a>