

MILLIN CAPSULE FORCEPS

angled; serrated

24.5cm (9.8")




SKU: 170-073-23

Categories: [Open Surgery](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Capsule Forceps, Ligature Carrying Forceps
Aesculap	EF166
KLS Martin	31-461-23
AUTHORS	MILLIN
Medicon	51.97.01



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Capsule Forceps, Ligature Carrying Forceps
Aesculap	EF166
KLS Martin	31-461-23
AUTHORS	MILLIN

Medicon

51.97.01