

## MILLIN CAPSULE FORCEPS

angled; serrated

24.5cm (9.8")




**SKU:** 170-073-23

**Categories:** [Open Surgery](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Capsule Forceps, Ligature Carrying Forceps</a>
<b>Aesculap</b>	EF166
<b>KLS Martin</b>	31-461-23
<b>AUTHORS</b>	<a href="#">MILLIN</a>
<b>Medicon</b>	51.97.01



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Capsule Forceps, Ligature Carrying Forceps</a>
<b>Aesculap</b>	EF166
<b>KLS Martin</b>	31-461-23
<b>AUTHORS</b>	<a href="#">MILLIN</a>

**Medicon**

51.97.01