

MINI-HOHMANN BONE LEVER

6mm

16.0cm (6.4")


SKU: 190-571-06

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Bone Levers |
| Aesculap | FK117 |
| KLS Martin | 23-539-06 |
| AUTHORS | MINI-TOWNSEND |
| Medicon | 78.35.11 |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Bone Levers |
| Aesculap | FK117 |
| KLS Martin | 23-539-06 |
| AUTHORS | MINI-TOWNSEND |

Medicon

78.35.11