

MIXTER HAEMOSTATIC FORCEPS

Fig. 3; curved; serrated

28.0cm (11.2")




SKU: 040-479-28

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| KLS Martin | 13-495-28 |
| AUTHORS | MIXTER |
| Medicon | 17.12.27 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| KLS Martin | 13-495-28 |
| AUTHORS | MIXTER |
| Medicon | 17.12.27 |

SPECIFICATIONS

[curved](#), [serrated](#)