

MIXTER HAEMOSTATIC FORCEPS

curved; serrated

15.5cm (6.2")



SKU: 040-461-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| KLS Martin | 13-473-15 |
| AUTHORS | MIXTER |
| Medicon | 17.07.36 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| KLS Martin | 13-473-15 |
| AUTHORS | MIXTER |
| Medicon | 17.07.36 |

SPECIFICATIONS

[curved](#), [serrated](#)