

## MODELL USA DRESSING FORCEPS

straight; serrated; narrow

13.0cm (5.2")

**SKU:** 020-012-13

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-146-13
<b>AUTHORS</b>	<a href="#">MODELL USA</a>
<b>Medicon</b>	06.34.53
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a>
<b>KLS Martin</b>	12-146-13
<b>AUTHORS</b>	<a href="#">MODELL USA</a>
<b>Medicon</b>	06.34.53

## SPECIFICATIONS

[serrated](#), [straight](#)