

MODELL USA DRESSING FORCEPS

straight; serrated; narrow

15.5cm (6.2"); width of jaw 1.9 +/-0.2


SKU: 020-012-15

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| KLS Martin | 12-146-15 |
| AUTHORS | MODELL USA |
| Medicon | 06.34.55 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| KLS Martin | 12-146-15 |
| AUTHORS | MODELL USA |
| Medicon | 06.34.55 |

SPECIFICATIONS

[serrated](#), [straight](#)