

## MODELL USA DRESSING FORCEPS

straight; serrated; narrow

15.5cm (6.2"); width of jaw 1.9 +/-0.2

**SKU:** 020-012-15

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>  |
| <b>TYPE</b>           | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>KLS Martin</b>     | 12-146-15   |
| <b>AUTHORS</b>        | <a href="#">MODELL USA</a>  |
| <b>Medicon</b>        | 06.34.55  |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a>               |



### ZUSÄTZLICHE INFORMATIONEN

|                   |   |
|-------------------|---|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Dressing Forceps</a> , <a href="#">Tissue Forceps</a> |
| <b>KLS Martin</b> | 12-146-15   |
| <b>AUTHORS</b>    | <a href="#">MODELL USA</a>  |
| <b>Medicon</b>    | 06.34.55  |

## SPECIFICATIONS

[serrated](#), [straight](#)