

MODELL USA DRESSING FORCEPS

straight; serrated

13.0cm (5.2"); width of jaw 3.0 +/-0.2

SKU: 020-006-13

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-13
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-13
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight