

MODELL USA DRESSING FORCEPS

straight; serrated

13.0cm (5.2"); width of jaw 3.0 +/-0.2

SKU: 020-006-13

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| KLS Martin | 12-140-13 |
| AUTHORS | MODELL USA |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| KLS Martin | 12-140-13 |
| AUTHORS | MODELL USA |
| SPECIFICATIONS | serrated , straight |