

MODELL USA DRESSING FORCEPS

straight; serrated

20.0cm (8"); width of jaw 3.9 +/-0.2

SKU: 020-006-20

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-21
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-21
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight