

MODELL USA DRESSING FORCEPS

straight; serrated

25.0cm (10"); width of jaw 4.0 +/-0.2

SKU: 020-006-25

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-25
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-140-25
AUTHORS	MODELL USA
SPECIFICATIONS	serrated , straight