

MODELL USA DRESSING FORCEPS

straight; serrated

25.0cm (10"); width of jaw 4.0 +/-0.2

SKU: 020-006-25

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| KLS Martin | 12-140-25 |
| AUTHORS | MODELL USA |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| KLS Martin | 12-140-25 |
| AUTHORS | MODELL USA |
| SPECIFICATIONS | serrated , straight |