

MOLLISON SELFRETAINING RETRACTOR

sharp; 4X4 teeth

15.0cm (6")




SKU: 070-674-15

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV088
KLS Martin	15-754-15
AUTHORS	MOLLISON
Medicon	24.26.30
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV088
KLS Martin	15-754-15

AUTHORS

[MOLLISON](#)

Medicon

24.26.30

SPECIFICATIONS

[sharp, toothed](#)