


## MONOPOLAR CABLE FOR

Martin units; 5m; Umax: 3 KVp

instrument and unit side Ø 4.0mm



**SKU:** 287-815-50

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Monopolar Cables</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Monopolar Cables</a>