


MONOPOLAR CABLE FOR ELECTRODES AND LOOPS

Ø 8mm; 3.0m; to Valleylab/conmed

equipment



SKU: 582-532-03

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Monopolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Monopolar Cables