


## MONOPOLAR CABLE FOR ELECTRODES AND LOOPS

Ø 4mm; 3.0m for equipm. with banana plug



**SKU:** 582-533-04

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Monopolar Cables</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Cables</a>
<b>TYPE</b>	<a href="#">Monopolar Cables</a>