


MONOPOLAR CABLE FOR ELECTRODES AND LOOPS

Ø 4mm; 3.0m for equipm. with banana plug



SKU: 582-533-04

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Monopolar Cables

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Cables
TYPE	Monopolar Cables