

MOUSE TEETH FORCEPS

opening to both sides; Ø1.2mm; FL 600mm

4 Fr.; rigid; with LuerLock; with handle



SKU: 583-480-12LL

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--|
| UNIT | pcs. |
| RANGE | Uretero-Renoscropy |
| TYPE | Semi-rigid Instruments |



ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--|
| UNIT | pcs. |
| RANGE | Uretero-Renoscropy |
| TYPE | Semi-rigid Instruments |