

MOUSE TEETH FORCEPS

opening to both sides; Ø1.2mm; FL 600mm

4 Fr.; rigid; with LuerLock; insert only



with Luer Lock

SKU: 583-480-12LLA

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Uretero-Renoscapy](#)

TYPE [Semi-rigid Instruments](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Uretero-Renoscapy](#)

TYPE [Semi-rigid Instruments](#)