

MOUSE TEETH/IUD FORCEPS

FL 600 mm; semi-rigid; 5 Fr.

double action; with LuerLock;insert only



with Luer Lock

SKU: 583-480-16LLA

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Uretero-Renoscapy](#)

TYPE [Semi-rigid Instruments](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Uretero-Renoscapy](#)

TYPE [Semi-rigid Instruments](#)