

## MOUSE TEETH/IUD FORCEPS

FL 600 mm; semi-rigid; 5 Fr.

double action; with LuerLock;insert only



with Luer Lock

**SKU:** 583-480-16LLA

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uretero-Renoscapy](#)

**TYPE** [Semi-rigid Instruments](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uretero-Renoscapy](#)

**TYPE** [Semi-rigid Instruments](#)