


MÜLLER APPROXIMATOR

Ø 2.0mm; closing force 25-30g



SKU: 026-530-02

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Clips, Approximators |
| AUTHORS | MUELLER |
| Medicon | 56.78.02 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Clips, Approximators |
| AUTHORS | MUELLER |
| Medicon | 56.78.02 |