

NASAL-BILE-DUCT DRAINAGE-SET, TYPE STAND

with 0.035" guide wire, 5fr. nasal

tube connector, funnel 3-way stopcock



SKU: PK-KNBS5

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Single Use Instruments](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Single Use Instruments](#)