

## NEEDLE CASE

perforated

65x42x8.0mm



**SKU:** 250-304-65

**Categories:** [Open Surgery](#), [Various Items](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>KLS Martin</b>	56-302-65

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>KLS Martin</b>	56-302-65