


## NEEDLE CASE

perforated

210x125x17mm



**SKU:** 250-310-20

**Categories:** [Open Surgery](#), [Various Items](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>Aesculap</b>	BL943
<b>KLS Martin</b>	55-162-20

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>Aesculap</b>	BL943
<b>KLS Martin</b>	55-162-20