

## NEEDLE CASE

perforated

50x25x4.5mm



**SKU:** 250-302-50

**Categories:** [Open Surgery](#), [Various Items](#)

## ZUSÄTZLICHE INFORMATIONEN

|                   |                       |
|-------------------|-----------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Bowls</a> |
| <b>Aesculap</b>   | BL910                 |
| <b>KLS Martin</b> | 56-302-50             |
| <b>Medicon</b>    | 46.37.46              |

## ZUSÄTZLICHE INFORMATIONEN

|                   |                       |
|-------------------|-----------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Bowls</a> |
| <b>Aesculap</b>   | BL910                 |
| <b>KLS Martin</b> | 56-302-50             |
| <b>Medicon</b>    | 46.37.46              |