


## NEEDLE CASES COMPLETE

with interior box

Ø 65.0mm



**SKU:** 250-300-00

**Categories:** [Open Surgery](#), [Various Items](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>Aesculap</b>	BL930
<b>KLS Martin</b>	55-309-65
<b>Medicon</b>	46.37.50

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>Aesculap</b>	BL930
<b>KLS Martin</b>	55-309-65
<b>Medicon</b>	46.37.50