

NEEDLE HOLDER Ø 5MM, CONVEX/ CONCAVE

with Aluminium handle, teflon coated

WL 450mm, LL-Connection



SKU: 614-120-54

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#)

TYPE [LOFTLINE](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Needle Holders](#)

TYPE [LOFTLINE](#)