

## NEEDLE HOLDER Ø 5MM, CONVEX/ CONCAVE

with Aluminium handle, teflon coated

WL 330mm, LL-Connection



**SKU:** 614-120-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Needle Holders](#)  
**TYPE** [LOFTLINE](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Needle Holders](#)  
**TYPE** [LOFTLINE](#)