

NEEDLE HOLDER Ø 5MM, CONVEX/ CONCAVE

with Aluminium handle, teflon coated

WL 330mm, LL-Connection



SKU: 614-120-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Needle Holders](#)
TYPE [LOFTLINE](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Needle Holders](#)
TYPE [LOFTLINE](#)