

NEIVERT NEEDLE HOLDER

straight; smooth; special groove

13.0cm (5.2")




SKU: 050-046-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
KLS Martin	20-338-13
AUTHORS	NEIVERT
Medicon	10.08.12
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
KLS Martin	20-338-13
AUTHORS	NEIVERT
Medicon	10.08.12

SPECIFICATIONS

[straight](#)