

## NEPHROSCOPE; 45°DOUBLE ANGLED COMPACT

18.8 Charr. oval (6.9mm x 5.6mm)

view 15°, working lenght 220mm

ACMI- Wolf- and Storz adapters


**SKU:** 601-135-45

**Categories:** [Endoscopy](#), [Equipment](#), [Urology](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">Telescopes</a> |
| <b>SCOPE</b> | <a href="#">Urology</a>    |



### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">Telescopes</a> |
| <b>SCOPE</b> | <a href="#">Urology</a>    |