

NEPHROSCOPY SHEATH; 23 FR.

FL 220mm; for Compact-Nephroscopes

601-135-45 + 601-135-90; 15°

colour code: white



SKU: 630-423-12

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Nephroscopy
TYPE	Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Nephroscopy
TYPE	Sheaths