

NITINOL-BASKET , TIPLESS ROUND

hexagonal; 25 X 40mm; 4 wires

2.6mm; 230cm



SKU: PK-1NRPK2625230K-Z

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)