

OBTURATOR VISUAL FOR RESECTOSCOPY

24 / 26Fr.; with 5Fr.working channel

for 4mm/30° telescope

SKU: 630-230-24

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators